


# TSLog 21



Deepwave system management software for continuous monitoring of environments, warehouses, cold rooms, fridge, ultra freezers, laboratories. It allows you to view all the data received, print the report, export it to Excel, manage remote alarms and users enabled to view and receive notifications.

Compatible with HACCP, 21 CFR Part11, Annex 11, GAMP 5 and validated in the pharmaceutical, cosmetic, health, food sectors.



## Main features

- Multi-level and multi-user
- Access through Login and Password
- Visualization of the data in graphical and table format
- Printable data in customizable reports and exportable to excel
- Management of remote alarms for multiple receivers and with indication of availability
- Compatible with 21 CFR PART 11, Annex 11, GAMP 5, HACCP
- Validable in healthcare, pharmaceutical and cosmetic fields
- Audit trail
- Electronic signature

## Plus

- Simple to configure: with a single software you manage the sensors, data loggers, data and alarms
- Easy calibration with Smart Sensors
- Management of availability shifts for alarms
- Scalable to multiple users: decide who can see what
- Environments map viewable on full screen with indication of read values and out of range values

## The system

The system is composed of:

- Lo.Ra. radio modules (Long Range) with Smart Sensor
- Gateway Lo.Ra Helios receiver
- Software 21 CFR Part 11 TSLog 21

## Accessories

- Helios - LoRa Gateway

## Applications



Pharmaceutical



Cosmetics



Healthcare



Medical



Warehouses



Cold Chain



Environment



Food & Beverages



Laboratories

